foam silicone

INSTRUCTIONS FOR USE

Kliniderm® Foam Silicone range consists of the following product variants:

- Kliniderm® Foam Silicone
- Kliniderm® Foam Silicone Border
- Kliniderm® Foam Silicone Lite
- Kliniderm® Foam Silicone Lite Border
- Kliniderm® Foam Silicone Heel
- Kliniderm® Foam Silicone Sacrum

Not made with natural rubber latex or DEHP. Medical device.

Indications

Kliniderm® Foam Silicone dressings are indicated for moderate to heavily exuding chronic and acute wounds and may be used throughout the healing process on the following wounds: pressure ulcers, venous leg ulcers, diabetic foot ulcers, lacerations, abrasions, skin tears, donor sites, postoperative surgical wounds, superficial and partial thickness burns.
Kliniderm® Foam Silicone dressings are suitable for use under compression bandaging.
Kliniderm® Foam Silicone Heel and Sacrum variant are ideal for use on the heel, sacrum and similar anatomical areas.

Contra-indications

None known.

Directions for use

Application:

- Select a suitable size of Kliniderm® Foam Silicone so that the dressing overlaps the wound margins by approximately 2cm/0.8". Non-bordered dressings can be cut to size if necessary.
- Remove the release film before use and apply the adherent side to the wound.
- If necessary, Kliniderm® Foam Silicone can be fixated with a bandage or other fixation.

IFU_KWC041_07-01

11-05-2023 / 0_040102

0123

Zhejiang Longterm Medical Technology Co., Ltd.

ShuangShan Road 277
Fuxi Street, Deqing County
313200 Huzhou City, Zhejiang Province
People's Republic of China
www.longterm-health.com

EU Authorized Rep.
Lotus NL BV
Koningin Julianaplein 10
1e verd, 2595 AA
The Hague, The Netherlands

- Do not use on patients with a known sensitivity to silicone, polyurethane films, foams or acrylic adhesives.
- Do not use on patients with a known sensitivity to silicone, polyurethane films, foams or acrylic adhesives.
- Do not re-use. Re-use can cause serious infections.
- Do not re-sterilize.
- For external use only.
- For use under the guidance of a healthcare professional.

In case a serious incident has occurred in relation to the device: Inform Medeco BV as listed on the packaging, your local distributor and your local healthcare authority.

Mediq Suisse AG
Rosengartenstrasse 25
CH - 8806 Bubikon
Switzerland

Medeco B.V.
Brandpuntlaan Zuid 14
2685 NZ Bleiswijk
The Netherlands

www.medeco.org

ENGLISH

Product Description
Kliniderm® Foam Silicone dressings are sterile wound dressings made from polyurethane foam. The wound contact layer is coated with a layer of soft silicone that does not stick to the wound surface and reduces trauma to delicate new tissue upon removal. This soft silicone layer is also slightly tacky, which facilitates application and retention of the dressing to intact skin, but does not cause epidermal stripping or pain upon removal. The center of the dressing is an absorbent polyurethane foam which can lock wound exudate into the product, preventing the drying out of the wound bed and maceration of the surrounding skin. The outer layer of the foam has a high permeability to allow effective vapor transfer, is waterproof and provides a bacterial barrier. As Kliniderm® Foam Silicone maintains a moist wound environment, supporting debridement, there might be an initial increase in the wound size. This is normal and to be expected.

Kliniderm® Foam Silicone range consists of the following product variants:

- Kliniderm® Foam Silicone
- Kliniderm® Foam Silicone Border
- Kliniderm® Foam Silicone Lite
- Kliniderm® Foam Silicone Lite Border
- Kliniderm® Foam Silicone Heel
- Kliniderm® Foam Silicone Sacrum

Not made with natural rubber latex or DEHP. Medical device.

Indications

Kliniderm® Foam Silicone dressings are indicated for moderate to heavily exuding chronic and acute wounds and may be used throughout the healing process on the following wounds: pressure ulcers, venous leg ulcers, diabetic foot ulcers, lacerations, abrasions, skin tears, donor sites, postoperative surgical wounds, superficial and partial thickness burns.
Kliniderm® Foam Silicone dressings are suitable for use under compression bandaging.
Kliniderm® Foam Silicone Heel and Sacrum variant are ideal for use on the heel, sacrum and similar anatomical areas.

Contra-indications

None known.

Directions for use

Application:

- Select a suitable size of Kliniderm® Foam Silicone so that the dressing overlaps the wound margins by approximately 2cm/0.8". Non-bordered dressings can be cut to size if necessary.
- Remove the release film before use and apply the adherent side to the wound.
- If necessary, Kliniderm® Foam Silicone can be fixated with a bandage or other fixation.

Dressing change and removal:

- Kliniderm® Foam Silicone dressings can remain in place for up to 7 days, dependent on the condition of the patient and the level of exudate.
- Replace Kliniderm® Foam Silicone when the dressing has reached its absorbent capacity or whenever good wound practice dictates that the dressing should be changed.
- Gently remove the dressing from the wound bed and skin (by supporting the surrounding skin) and dispose according to local procedures and guidelines. If removal appears to be difficult, saturate the dressing with sterile saline solution.
- If required, irrigate the wound site according to good wound care practice.

Precautions

- In case of signs of infection, consult a healthcare professional for adequate infection treatment.
- Do not use with oxidising solutions such as hypochlorite or hydrogen peroxide, as these can break down the absorbent component of the dressing.
- Do not use on patients with a known sensitivity to silicone, polyurethane films, foams or acrylic adhesives.
- Sterile if pouch and seal are intact. Do not use if damaged.
- Do not re-use. Re-use can cause serious infections.
- Do not re-sterilize.
- For external use only.
- For use under the guidance of a healthcare professional.

In case a serious incident has occurred in relation to the device: Inform Medeco BV as listed on the packaging, your local distributor and your local healthcare authority.

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