
klinidern®
super fiber

CMC Gelling Fiber
Dressing Standard
INSTRUCTIONS FOR USE

ENGLISH

Product Description

Klinidern® Super Fiber dressings are single use, sterile nonwoven wound dressings, composed of 100% sodium carboxymethylcellulose fibers. This conformable and highly absorbent dressing absorbs wound exudate and transforms into a cohesive transparent gel, maintaining a moist environment to support the wound healing process and minimizing pain and discomfort to the patient during dressing changes. When the dressing is in situ, it absorbs and retains exudate to lock away harmful components contained within the wound exudate, such as bacteria and proteinases, reducing the risk of wound infection and maceration.

Not made with natural rubber latex or DEHP.

Indications

Klinidern® Super Fiber dressings are indicated for the management of moderate to heavily exuding chronic and acute wounds and ulcers of the Epidermis or Dermis dermis. The dressings may be used through the healing process on leg ulcers, pressure ulcers (stage II-IV), diabetic ulcers, surgical wounds (e.g., post-operative, wounds left to heal by secondary or tertiary intention), partial thickness burns and traumatic wounds (e.g., abrasions and lacerations).

This device is for short-term use (single use times2days and accumulated use times28 days).

Contra-Indications

- Do not use on surgical implantation.
- Do not use as a treatment for arterial bleeding or on heavy bleeding wounds.
- Do not use on patients who are hypersensitive to, or have had an allergic reaction to the dressing or its components.

Adverse reaction

Very occasional, potential risks associated with the use of Klinidern® Super Fiber dressings:

- Adverse reactions (blisters, hyper granulation, eczema, etc.)
- Wound infections
- Risk of maceration (fluid leakage)

Directions for use

- Application:
- Cleanse the wound in accordance with local procedures and guidelines. Ensure the peri-wound skin is dry.
 - Select a suitable size of Klinidern® Super Fiber so that the dressing overlaps the wound margins by at least 1 cm. The dressing may be cut or folded. When using the ribbon dressing, loosely pack the dressing to about 80% capacity leaving at least 2.5cm outside the wound for easy retrieval.
 - Gently apply the dressing to the wound and cover with an appropriate, moisture retentive secondary dressing, such as, but not limited to, transparent film dressing, hydrocolloid dressing, foam dressing or bandage. The secondary dressing should be a moisture retentive dressing, and the size of the dressing should be able to completely cover the wound area and overlap at least 2-3cm onto the skin surrounding the wound.

Dressing change and removal

- Klinidern® Super Fiber can remain in place for up to 7 days, dependent on the condition of the patient and the level of exudate. Cumulative use should not exceed 28 days.
- Replace Klinidern® Super Fiber when the dressing has reached its absorbent capacity or whenever good wound practice dictates that the dressing should be changed.
- Gently lift off the secondary dressing and remove the dressing from the wound bed and dispose according to local procedures and guidelines. If the dressing appears dry or when removal appears to be difficult, fully saturate the dressing with sterile saline solution.
- If required, irrigate the wound site with sterile saline solution according to good wound care practice for removal of gel residue.

Precautions

- In case of signs of infection, consult a healthcare professional for adequate infection treatment.
- Do not use on patients with a known sensitivity to the dressing or its components.
- Check if pouch and seal are intact; do not use if damaged.
- Do not re-use. Re-use can cause serious infections.
- Do not re-sterilise.
- For external use only.
- For use under the guidance of a healthcare professional.

Storage

The products should be stored under ambient conditions, in the original packaging. In a cool, dry and well-ventilated place, non-corrosive gas environment, away from direct sunlight, fire, and flammable materials.

In case a serious incident has occurred in relation to the device, inform Winner Medical Co., Ltd. as well as on the packaging, your local distributor and the competent healthcare authority.

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