kliniderm[®] foam phmb

INSTRUCTIONS FOR USE

ENGLISH

Product Description

Kliniderm[®] Non-Adhesive PHMB Foam dressing is a sterile wound dressing, consisting of a hydrophilic polyurethane foam with high absorption- and exudate retention capacity. The dressing is externally covered by a beige colored thin superficial film - acting as semi-permeable barrier - which optimizes the breathability, avoids the fluid leaking and protects the wound from external contamination and from the water. The wound contact layer side of the dressing has an particular open cell structure, with vertical absorption. This limits the side exudates dispersion and reduces the risk of wound maceration. The dressing is soft, elastic and extremely conformable, which ensures painless replacement of the dressing, even after some days of treatment. Kliniderm[®] Non-Adhesive PHMB Foam dressing contains 0.5% of the antimicrobial substance polyhexamethylene biguanide (PHMB) and 0,1 % B-Panthenol. The presence of PHMB into the pad ensures the microbiological control of the wound dressing, reducing the risks of infection and promoting the healing process. B-Panthenol improves the moisturizing of reepithelization tissues.

Indications

Kliniderm[®] Non-Adhesive PHMB Foam dressing is indicated for the treatment of chronic or acute wounds, moderate to heavily exuding. It can be applied on heavily exuding wound with ongoing infection under medical surveillance. The dressings may be used during the healing process on pressure ulcers, leg and foot ulcers, diabetic ulcers, and surgical wounds. Kliniderm[®] Non-Adhesive PHMB Foam is suitable for use under compression bandaging.

Contraindications

Kliniderm[®] Non-Adhesive PHMB Foam dressings are not indicated for use on the following: Dry or lightly exuding wounds, Individuals with a known sensitivity to polyurethane films, foams or PHMB.

Directions for use

Site Preparation
Dry the surrounding areas of the wound thoroughly before dressing application. Clean, disinfect and treat the wound according to the local protocol; apply the dressing on the wound and fixate with appropriate fixation devices.

Dressing Application
1. Select a size of Kliniderm[®] Non-Adhesive PHMB Foam that is slightly larger than the wound. If necessary the dressing can be cut (using sterile scissors) to the wound size.
2. Centre the dressing on the wound, ensuring the beige side is facing away from the wound and apply it gently to the wound site.

3. Fix the dressing in place with a suitable tape or fixation bandage to ensure the dressing remains in position. Ensure the tape is applied to the edges only and does not cover the dressing as this may affect its performance. Alternatively, if being used under compression bandaging, tape is not required as the bandaging will hold the dressing in situ.
4. Discard any remaining dressing material due to the risk of contamination.

Dressing Change and Removal
Kliniderm[®] Non-Adhesive PHMB Foam can remain in situ up to 7 days, dependent on patient condition and the level of exudate. Initially, it may be necessary to change the dressing every 24 hours or whenever good wound care practice dictates.
2. Gently remove the dressing from the wound bed and dispose of according to local procedures and guidelines.
3. If required, cleanse the wound site with a suitable wound cleanser prior to application of a new dressing.

Precautions

Do not use with oxidising solutions such as hypochlorite or hydrogen peroxide, as these can break down the absorbent polyurethane component of the dressing and can also break down the PHMB. In the event of clinical infection, PHMB does not replace the need for systemic therapy or other adequate infection treatment, consult the relevant healthcare professional. Suspend the application in case of worsening of the lesion and phenomena such as redness of the peri-wound tissue and/or irritation. Pediatric use and use on pregnant and lactating women requires strict medical supervision. Foam PHMB must not be used on newborn patients. Do not re-use in whole or in part, as it may compromise the ability and/or the performance of the dressing. Kliniderm[®] Non-Adhesive PHMB Foam dressings must not be ingested and must be kept away from children and animals.

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