kliniderm® foam silicone border flexible

INSTRUCTIONS FOR USE

KLINION. WOUND CARE

IFU_KWC063_01

02-09-2024 / 0_040101

CE 0123

Zhejiang Longterm Medical Technology Co., Ltd.
ShuangShan Road 277
Fuxi Street, Deying County
313200 Huzhou City, Zhejiang Province
People's Republic of China
www.longterm-health.com

EC REP
EU Authorized Rep.
Lotus NL BV
Koninklijke Julianaplein 10
1e verd., 2595 AA
The Hague, The Netherlands

EU Imported + Distributed by
Medeco B.V.
Brandpuntlaan Zuid 14
2665 NZ Bleiswijk
The Netherlands

Swiss Importer
Mediq Suisse AG
Rosengartenstrasse 25
CH - 8606 Bubikon
Switzerland

ENGLISH

Product Description

Kliniderm® Foam Silicone dressings are sterile wound dressings made from polyurethane foam. The wound contact layer is coated with a layer of soft silicone that does not stick to the wound surface and reduces trauma to delicate new tissue upon removal. This soft silicone layer is also slightly tacky, which facilitates application and retention of the dressing to intact skin, but does not cause epidermal stripping or pain upon removal. The center of the dressing is an absorbent polyurethane foam which can lock wound exudate into the product, preventing the drying out of the wound bed and maceration of the surrounding skin. The outer layer of the foam has a high permeability to allow effective vapor transfer, is waterproof and provides a bacterial barrier. As Kliniderm® Foam Silicone maintains a moist wound environment, supporting debridement, there might be an initial increase in the wound size. This is normal and to be expected.

Not made with natural rubber latex or DEHP. Medical device.

Indications

Kliniderm® Foam Silicone dressings are indicated for moderate to heavily exuding chronic and acute wounds and may be used throughout the healing process on the following wounds: pressure ulcers, venous leg ulcers, diabetic foot ulcers, lacerations, abrasions, skin tears, donor sites, postoperative surgical wounds, superficial and partial thickness burns. Kliniderm® Foam Silicone dressings are suitable for use under compression bandaging. Kliniderm® Foam Silicone Heel and Sacrum variant are ideal for use on the heel, sacrum and similar anatomical areas.

Contra-indications

None known.

Directions for use

Application:

1. Select a suitable size of Kliniderm® Foam Silicone so that the dressing overlaps the wound margins by approximately 2cm/0.8". Non-bordered dressings can be cut to size if necessary.
2. Remove the release film before use and apply the adherent side to the wound.
3. If necessary, Kliniderm® Foam Silicone can be fixed with a bandage or other fixation.

Dressing change and removal:

1. Kliniderm® Foam Silicone dressings can remain in place for up to 7 days, dependent on the condition of the patient and the level of exudate.
2. Replace Kliniderm® Foam Silicone when the dressing has reached its absorbent capacity or whenever good wound practice dictates that the dressing should be changed.
3. Gently remove the dressing from the wound bed and skin (by supporting the surrounding skin) and dispose according to local procedures and guidelines. If removal appears to be difficult, saturate the dressing with sterile saline solution.
4. If required, irrigate the wound site according to good wound care practice.

Precautions

- In case of signs of infection, consult a healthcare professional for adequate infection treatment.
- Do not use with oxidising solutions such as hypochlorite or hydrogen peroxide, as these can break down the absorbent component of the dressing.
- Do not use on patients with a known sensitivity to silicone, polyurethane films, foams or acrylic adhesives.
- Sterile if pouch and seal are intact. Do not use if damaged.
- Do not re-use. Re-use can cause serious infections.
- Do not re-sterilize.
- For external use only.
- For use under the guidance of a healthcare professional.

In case a serious incident has occurred in relation to the device: Inform Medeco BV as listed on the packaging, your local distributor and the competent authority.

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