kliniderm[®] foam
INSTRUCTIONS FOR USE

KLINION[®]

Medeco B.V.
Alexander Flemingstraat 2
NL-3261 MA Oud-Beijerland
The Netherlands

www.medeco.org

STERILE R

LATEX DEHP

MD

IFU_KWC056_06
02-10-2019 / 0_040101

ENGLISH

Product Description

Kliniderm Foam dressings are sterile wound dressings. The wound contact layer consists of hydrophilic, absorbent polyurethane foam. The outer layer is pink, smooth, waterproof and made of polyurethane film. It has a high permeability to allow effective vapour transfer and provides a bacterial barrier. Kliniderm Foam dressings maintain a moist wound environment to assist in the healing process of moderate to heavily exuding wounds. Kliniderm Foam dressings are available in three variants: adhesive/border, non-adhesive and heel.

Not made with natural rubber latex of DEHP. Medical device.

Indications

Kliniderm Foam dressings are indicated for moderate to heavily exuding chronic and acute wounds and may be used throughout the healing process on the following wounds: pressure ulcers, venous leg ulcers, diabetic foot ulcers, lacerations, abrasions, post-operative surgical wounds, superficial and partial thickness burns. Kliniderm Foam dressings are suitable for use under compression bandaging and may be used on cavity wounds as a secondary dressing. Kliniderm Foam Heel variant is ideal for use on the heel and similar anatomical areas.

Contra-indications

None known.

Directions for use

Application:

- Select a suitable size of Kliniderm Foam so that the dressing overlaps the wound margins by approximately 2 cm/0.8". Cut the dressing to size if necessary.
- Apply to the wound directly.
- In case of using the non-adhesive or heel variant of Kliniderm Foam, ensure the dressing is applied to the wound with the pink side up.
- In case of using a non-adhesive or heel variant of Kliniderm Foam, secure the dressing in place with an appropriate retention bandage or tape. The tape should only be applied to the edge of the dressing.

Dressing change and removal:

- Kliniderm Foam dressings can remain in

place for up to 7 days, dependent on the condition of the patient and the level of exudate.

- Replace Kliniderm Foam when the dressing has reached its absorbent capacity or whenever good wound practice dictates the dressing should be changed.
- Gently remove the dressing from the wound bed and dispose according to local procedures and guidelines. If removal appears to be difficult, saturate the dressing with sterile saline solution.
- If required, irrigate the wound site according to good wound care practice.

Precautions

- In case of signs of infection, consult a healthcare professional for adequate infection treatment.
- Do not use with oxidising solutions such as hypochlorite or hydrogen peroxide, as these can break down the absorbent component of the dressing.
- Do not use on patients with a known sensitivity to polyurethane films, foams or acrylic adhesives (border).
- Sterile if pouch and seal are intact. Do not use if damaged.
- Do not re-use. Re-use can cause serious infections.
- Do not re-sterilize.
- For external use only.
- For use under the guidance of a healthcare professional.

In case a serious incident has occurred in relation to the device: Inform Medeco BV as listed on the packaging, your local distributor and your local healthcare authority.

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