klinitulle® ointment compress

INSTRUCTIONS FOR USE

KLINION®

Medeco B.V.
Alexander Flemingstraat 2
NL-3261 MA Oud-Beijerland
The Netherlands

www.medeco.org

MD

ENGLISH

Product Description

Klinitulle® ointment compress is a nonwoven pad impregnated with a white paraffin ointment. The dressing is intended to be used as a barrier between the broken skin and the surrounding environment and will facilitate the wound healing by maintaining a moist wound environment. The ointment prevents adhering of the dressing to the wound. The dressing smoothes and protects newly formed granulation tissue and allows exudate to pass into a secondary absorbent dressing.

MD

Medical device.

Indications

Klinitulle® ointment compress is intended to be used in conjunction with a secondary absorbent dressing for superficial wounds and burns, such as minor burns and scalds, abrasions and lacerations, donor recipient graft sites, nail extractions and second-degree burns with an intact blister.

Contra-indications

- Do not use on heavily exuding wounds.

Directions for use

Application:

1. Make sure the wound area is clean and dry before using the dressing.
2. Select a suitable size of Klinitulle® ointment compress.
3. Remove the release film on both sides of the dressing before use.
4. Gently apply the dressing to the wound.
5. Select and cover the Klinitulle® ointment compress with a suitable secondary absorbent dressing.
6. If required, fixate the Klinitulle® ointment dressing and the secondary absorbent dressing with an appropriate bandage or tape.

Dressing change and removal:

1. Klinitulle® ointment compress can remain in place for up to 1 day, dependent on the condition of the patient, the clinical condition of the wound and the level of exudate. The longest cumulative use should be less than 30 days.
2. Replace Klinitulle® ointment compress at least once a day or whenever good wound practice dictates that the dressing should be changed.
3. Gently remove the dressing from the wound bed and dispose according to local procedures and guidelines. If removal appears to be difficult, saturate the dressing with sterile saline solution.
4. If required, irrigate the wound site according to good wound care practice.

Precautions

- In case of signs of infection, consult a healthcare professional for adequate infection treatment.
- Do not use on patients with a known sensitivity to the product itself or to its components.
- Sterile if pouch and seal are intact. Do not use if damaged.
- Do not re-use. Re-use can cause serious infections.
- Do not re-sterilize.
- For external use only.
- Paraffin-based products should be kept away from open flame as they are highly flammable.

In case a serious incident has occurred in relation to the device: Inform Medeco BV as listed on the packaging, your local distributor and your local healthcare authority.

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