


Kliniderm®
**foam silicone
border flexible**



Gentle • Flexible • Conformable

Kliniderm® foam silicone border flexible

Kliniderm® covers the advanced products of Klinion - a well - known wound care brand. For over 30 years, Klinion offers patients global good quality products for competitive prices. With a wide range, Kliniderm® product ensure optimal user comfort.

We are excited to announce the extension of our Kliniderm® Foam Silicone Border range with a new extra flexible product group: Kliniderm® Foam Silicone Border Flexible! This new product group consists of 7 sizes and offers extra flexibility and extensibility. The Kliniderm® Foam Silicone Border Flexible has unique Z- shaped cuts that allow for a more flexible dressing. This dressing can be used when extra flexibility is needed. Due to the special design the dressing is able to stay in place better during regular wear.

The transparent design allows the healthcare professional to check in on the level and type of exudate without disturbing the wound. This might guide the professional in dressing changes, or whenever good wound practice dictates that the dressing should be changed.

Indications

Kliniderm® Foam Silicone border flexible is indicated for moderate to heavily exuding chronic and acute wounds and may be used throughout the healing process on the following wounds: pressure ulcers, venous leg ulcers, diabetic foot ulcers, lacerations, abrasions, skin tears, donor sites, postoperative surgical wounds, superficial and partial thickness burns.

Gentle, flexible and conformable



Extra flexibility

Unique Z- shaped cuts that allow for a more flexible dressing and greater extensibility



Transparent design

Allows for better visibility on the wound during the healing process



Absorbs and manages exudate

To protect surrounding skin and support moist wound healing



Gentle silicone adhesive

Helps minimise the risk of pain and damage to the wound at dressing changes




Skin friendly

Reliable fixation and safe for fragile skin, improving patient compliance



Up to 7 days wear time

To minimise dressing changes



Effective under compression

Absorbs under compression bandages and hosiery



Protective

Providing a barrier for microorganisms



A MEDIQ COMPANY

Medeco B.V.
Brandpuntlaan Zuid 14
2665 NZ Bleiswijk
The Netherlands



www.medeco.org

MED_01_20241120