

# kliniderm<sup>®</sup> **silicone wound contact layer**

---

**Gentle contact**

# Kliniderm®.

Kliniderm® covers the advanced products of Klinion - a well - known wound care brand.  
For over 30 years, Klinion offers patients worldwide **good quality products** for **competitive prices**.  
With a wide range, Kliniderm® products ensure optimal user comfort.

## Indications for use.

Kliniderm® silicone wound contact layer is indicated for chronic and acute wounds including; skin tears, traumatic wounds, cuts, abrasions and lacerations, blisters, superficial and partial thickness burns, surgical wounds, pressure ulcers, diabetic foot ulcers and leg ulcers.

*Kliniderm silicone wound contact layer  
provides significant savings  
Ask to find out more!*



## Gentle contact, ensuring patient comfort.

Kliniderm® silicone wound contact layer consists of a polyurethane film coated with a soft silicone layer. Its porous structure allows exudate to pass through easily into an outer absorbent dressing such as Kliniderm® superabsorbent.



### Gentle silicone adhesive

*helps minimise the risk of pain and damage to the wound at dressing changes*



### Skin friendly

*reliable fixation and safe for fragile skin improving patient compliance*



### Up to 14 days wear time

*allowing the wound to remain undisturbed when the secondary dressing is changed*




### One sided adhesive contact layer

*ensures easy handling and prevents secondary dressings from sticking*



### Transparency

*allows for wound inspection*



### Reduces risk of maceration

*perforations allow exudate to pass into secondary dressing*



### Conformable and comfortable

*for the patient, helping to improve compliance*



### Easy application and removal

*helps minimise wound disturbance*



### Can be cut to shape

*to suit clinical requirements*



A MEDIQ COMPANY

### Medeco B.V.

Alexander Flemingstraat 2, NL-3261 MA Oud-Beijerland  
The Netherlands



[www.medeco.org](http://www.medeco.org)